

WAM™ Arm

Quick Start Guide

This Quick Start Guide steps you through connecting to the WAM system, zero calibration, gravity calibration and running a simple program. This can be used in conjunction with the additional resources available at <http://support.barrett.com/wiki/WAM> for new users.

Lab Space Setup

Requirements for operation:

- 100-240 VAC, 50-60 Hz
- Stable, flat surface to mount WAM (such as the included WAM stand)
- An unobstructed workspace with a 1 meter radius around the WAM
- A notebook PC with a wireless card for logging in to the WAM PC

Install and Connect to the WAM

- Mount WAM on a suitable platform (see <http://support.barrett.com/wiki/WAM/HardwareSetup>)
- Arrange and connect the electrical cables as shown in Figure 1
- *Note:* If you have a Display Pendant, you may connect it to the Display Pendant port. Otherwise, plug the Pendant Jumper into that port.
- Mount the provided wireless access point (AP) to the WAM backplate and connect it to the network port of the WAM as shown in Figure 2
- Turn on the WAM's power supply
 - Safety system comes online
 - Safety system loads defaults from non-volatile RAM/EEPROM
 - Safety system mode is set to JOINT velocity limiting
 - Safety system clears pendants and exercises pendant lights
 - Safety system observes critical bus undervoltage, critical fault (normal)
 - Safety system resistively brakes the WAM's power bus
- Disconnect any Ethernet cables from your notebook PC
- Choose "WAMxxx" as your wireless access point SSID, where xxx is your WAM's serial number
- Windows users:
 - Download the free PuTTY SSH client
 - Download Barrett's wamdiscover.exe application http://web.barrett.com/support/WAM_Installer/wamdiscover.exe
 - Put the two applications next to each other on your desktop
- Linux users:
 - Download Barrett's wamdiscover.py Python script for Linux: http://web.barrett.com/support/WAM_Installer/wamdiscover.py
- Run wamdiscover and select your WAM to log in using ssh
 - Username: robot
 - Password: WAM




Figure 1




Figure 2



Barrett Technology™ Inc.

Calibrate the Robot

- Make sure both pendant E-STOP buttons (red) are reset (up)
- Press <Shift+Reset/Idle> on the Control Pendant
 - Safety system clears the voltage fault
 - Safety system releases the resistive braking
 - Safety system turns on the power bus
 - Motor controllers come online and wait for “Initialize” command
 - There is no WAM braking at this stage
 - Safety system enumerates motor controllers
 - *** Safety system is fully operational now ***
 - Safety system registers heartbeat error, displays warning
- Place the WAM into the home (folded) position (see Figure 3)
 - Contact support@barrett.com if your WAM is mounted in a different orientation
- Zero-Calibrate the WAM using the following command:

```
$ bt-wam-zerocal
```

 - Program sends “Initialize” command to motor controllers
 - Motor controllers load defaults from flash/EEPROM
 - Motor controllers apply FET braking at motors
 - Program sets the motor controllers’ initial positions
 - Safety system is set to CARTESIAN velocity limiting
 - Program sets the safety system velocity and torque limits
 - Program requests positions from the motor controllers
 - Motor controllers send positions to the program
 - Safety system calculates/limits elbow & arm end-tip velocities
 - Program sends zero-value torques to the motor controllers
 - Safety system records latest torques
 - Safety system clears heartbeat error (it sees periodic torques & positions)
 - *** Pendant should now display all OK ***
 - Press Shift-Activate on the Control Pendant
 - If there are no warnings or faults, safety system tells controllers to “Activate”
 - Motor controllers release their FET braking and apply received torques
 - *** Motor controllers are actively commutating zero torque now ***
 - Please follow the on-screen instructions. Additional detailed instructions can be found here:
 - http://web.barrett.com/support/WAM_Documentation/zero_calibration_procedure.pdf
- Gravity-Calibrate the WAM using the following command:

```
$ bt-wam-gravitycal
```

 - Please follow the on-screen instructions to complete the WAM’s gravity calibration




Figure 3

Build and Run the Demo Application

- ```
$ cd libbarrett_examples
```
- ```
$ cmake . (cmake <space> <period>)
```
- ```
$ make ex08_teach_and_play
```
- ```
$ ./ex08_teach_and_play
```

 - Host program calculates and applies gravity compensation torques for the WAM
 - WAM should be free-floating now
- Press <Enter> to begin teaching the WAM a trajectory
 - The application is recording the current joint positions at 500 Hz
- Grab the WAM and move it slowly through a trajectory of your choice
- Press <Enter> to finish teaching the trajectory
- Press <Enter> to play back the recorded trajectory
 - The host program commands the robot to follow the recorded joint trajectory

Congratulations! You have set up your WAM Arm!